

CuroCell® S.A.C.



Seat cushion for support & comfort

CuroCell® S.A.C.

CuroCell® S.A.C. is a static air seat cushion for comfort, prevention and as an aid in the treatment of pressure ulcers up to category 2⁽¹⁾.

Facts

Type: Static air seat cushion

Function: Prevention and treatment up to and including category 2⁽¹⁾

Recommended user weight: 0-160 kg

Material: Film, color gray

Sizes: 45x45x5 cm & 50x50x5 cm

Warranty: 1 year

Classification: ISO code 04 33 03 (pressure relieving)

Cleaning: Wipe with cleaning agent and/or disinfectants



It is pre-filled and easy to handle. The durable soft material for long term use provides a good comfort and support.

Features and Characteristics

- Air venting holes, allows transportation of heat and moisture
- Durable soft material for long term use
- Light and easy to use
- Low profile design
- Patient ready
- Pre-inflated
- Reduces tissue interface pressure

References

(1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell are trademarks of Care of Sweden AB.

© Care of Sweden AB, 2015. All rights reserved. Last update: 02/06/15.